



hsa-miR-4725-3p

5'	gugucucucuggagaccucgagccuuccaccaccaggaggagcuuccaugggcugugggaaggcgucagugucgggugaggggaacac	-3'	exp	
	((((((((.....(((.(((.((((((((((.....(((.(.....)).)).)).)).)).)).)).)).)).)).)).)).))	reads	mm	sample
.....	agaccucgagccuuccacccaU.....	1	1	hsa
.....	uggggaaggcgucagugucg.....	2	0	hsa
.....	uggggaaggcgucagugucgg.....	1	0	hsa
.....	uggggaaggcgucagugucggA.....	1	1	hsa
.....	uggggaaggcgucagugucggg.....	6	0	hsa
.....	uggggaaggcgucagugucgggA.....	5	1	hsa
.....	uggggaaggcgucagugucgggu.....	4	0	hsa