



hsa-miR-491-3p

5' -	agugacuuagcuggguuaguggggaacccuuccaugaggaguagaacacuc <u>cuuau</u> gcaagau <u>ucccuuac</u> cugggcuggguugg	-3'	exp	
(.((((((.((((.((((.((((.((((((((((.(...)))))))))).)))))))))).)))))).)))))).))	reads	mm	sample	
.....agugggga <u>U</u> ccuuccaugagg.....	1	1	hsa	
.....aguggggaacccuuccaugagg.....	3	0	hsa	
.....aguggggaacccuuccaugagga.....	13	0	hsa	
.....aguggggaacccuuccaugaggaA.....	1	1	hsa	
.....aguggggaacccuuccaugaggag.....	1	0	hsa	
.....aguggggaacccuuccaugaggaU.....	1	1	hsa	
.....aguggggaacccuuccaugaggaA.....	1	1	hsa	
.....aguggggaacccuuccaugaggaCu.....	1	1	hsa	
.....uggggaacccuuccaugaggagu.....	3	0	hsa	
..... <u>cuuau</u> gcaagau <u>ucccuuac</u>	2	0	hsa	